

Air-cooled diode-pumped Q-switched laser



Laser head Q2

Q2 series diode-pumped, fully air-cooled, Q-switched laser intended for wide range of applications that require high peak power pulses.

Our innovative water-free laser crystal end-pumping technology allows to produce Gaussian-like, low divergence laser beam. At the same time, Q2 is a versatile platform that can be configured in many ways.

When used with the **H-SMART stand-alone harmonic generator**, the system can provide up to **3 separate output ports** with different harmonics, extending the applicability of the Q2 lasers.



H-SMART harmonics generator with fiber coupler FC
Provides up to 5 harmonics, compatible with the Q2 lasers

Q2

FEATURES

- **1053 nm** or **1064 nm** wavelength from **Nd:YLF** or **Nd:YAG** laser crystals respectively
- Up to **80 mJ** pulse energy and up to **2 W** average power
- Up to **100 Hz** pulse repetition rate
- **Air-cooled** (water-free)
- **5–10 ns** pulse duration
- Optional **< 3 ns** pulse duration at up to **40 mJ** pulse energy (short cavity version)
- Weight **< 9 kg** incl. heatsink, controller & AC/DC adapter
- Guaranteed **> 2 G shot** lifetime of pump-diodes
- **Built-in sync pulse generator** for triggering of user equipment
- **Remote monitoring and control** via built-in **Ethernet** interface

OPTIONAL EQUIPMENT

Stand alone:

- H-SMART series stand-alone **2nd, 3rd, 4th or 5th** harmonic generators
- **Two-channel pulse generator** for triggering
- **Air-purging unit** for long lifetime of UV optics

Attachable:

- **2nd harmonic generator**, model SHG
- **Beam guiding** module
- **Motorized attenuator** for fundamental wavelength beam
- **Pulse energy monitor** with analog and/or digital output

APPLICATIONS

- Light Detection And Ranging (LIDAR) (Metrology, Military, Astronomy and etc.)
- Spectroscopy (Laser Induced Breakdown (LIBS), Light Induced Fluorescence (LIFS), Flash photolysis and etc.)
- Laser ablation (marking, LCD repair, trimming, scribing and etc.)
- Medical (Ophthalmology, Dermatology, Photoacoustic Imaging (PAI) and etc.)
- Pulsed Laser Deposition (PLD)

SPECIFICATIONS at 10 Hz pulse repetition rate ¹⁾

| Model | Q2 | | | | |
|--------------------------------------|-------------------------------------|----------|----------|----------|---------|
| | -B10 | -C10 | -D10 | -E10 | -F10 |
| Pulse repetition rate ²⁾ | 10 Hz | | | | |
| Wavelength, nm | 1064/1053 nm | | | | 1053 nm |
| Pulse energy ³⁾ | 8/10 mJ | 16/20 mJ | 32/40 mJ | 60/65 mJ | 80 mJ |
| Typical pulse duration ⁴⁾ | < 8 ns | | < 7 ns | | < 5 ns |
| Beam divergence ⁵⁾ | < 1.5 mrad | | < 1 mrad | | |
| Typical beam diameter ⁶⁾ | 1.5 mm | | 3.0 mm | | 4.0 mm |
| Pulse energy stability | < 0.5 % RMS | | | | |
| Power drift | ± 3.0 % | | | | |
| Beam profile | bell-shaped, > 80 % fit to Gaussian | | | | |
| Jitter | < 0.5 ns RMS | | | | |
| Polarization | linear, horizontal | | | | |

Optional harmonics generator ⁷⁾

| Pulse energy ³⁾ | -B10 | -C10 | -D10 | -E10 | -F10 |
|----------------------------|------------|----------|----------|----------|-------|
| 532/526.5 nm | 4/5 mJ | 8/10 mJ | 16/20 mJ | 30/32 mJ | 40 mJ |
| 355/351 nm | 2.4/3 mJ | 4.5/6 mJ | 9/12 mJ | 18/20 mJ | 24 mJ |
| 266/263 nm | 1.2/1.5 mJ | 2.4/3 mJ | 5/6 mJ | 9/10 mJ | 12 mJ |
| 213/211 nm | 0.5 mJ | 1 mJ | 2 mJ | 3.5 mJ | 5 mJ |


Optional attenuator ⁸⁾

| | |
|--------------------|------------|
| Transmission range | 0.5 – 95 % |
|--------------------|------------|

Operating requirements

| | |
|---------------------------|--|
| Cooling requirements | air-cooled (water-free) |
| Ambient temperature | 15 – 30 °C |
| Relative humidity | 10 – 80 % (non-condensing) |
| Mains voltage | 90 – 230 VAC, single phase, 47 – 63 Hz ⁹⁾ |
| Average power consumption | 30 W 40 W 50 W 60 W |

1. Due to continuous improvements all specifications are subject to change. Unless stated otherwise all specifications are measured at fundamental wavelength and maximum pulse repetition rate. The parameters marked typical are not specifications. They are indications of typical performance and will vary with each unit we manufacture.
2. Factory-set pulse repetition rate is fixed at max repetition rate shown in the table. In internal triggering mode repetition rate can be divided by integer number down to f/2, f/3, f/4... 1 Hz.
3. When pulse energy is presented in xx/yy format, first number is for 1064 nm wavelength version, second – for 1053 nm version.
4. At FWHM level at fundamental wavelength, measured with 350 ps rise time photodiode. Short pulse duration version is available, with pulse duration shorter by approx 50 %. Inquire for detailed specifications.
5. Full angle measured at the 4σ level.
6. Beam diameter is measured 20 cm from laser output at the 4σ level.
7. Q2 is compatible with our attachable second harmonic generator (model SHG) and all models of stand-alone H-SMART harmonics generator. Pulse energies presented here are maximum values. Please refer to harmonic generator datasheets for detailed specifications.
8. Motorized attenuator intended to be attached to the laser housing. Transmission can be changed remotely through laser web-server control interface.
9. Laser can be powered from an appropriate 12 or 28 VDC power source. Please inquire for details.



Laser head Q2 and H-SMART harmonics generator on BSP2 base plate

SPECIFICATIONS at 20 – 33 Hz pulse repetition rate ¹⁾

| Model | Q2 | | | | | | |
|--------------------------------------|-------------------------------------|----------|----------|----------|---------|--------------|---------|
| | -B20 | -C20 | -D20 | -E20 | -F20 | -D33 | -E33 |
| Pulse repetition rate ²⁾ | 20 Hz | | | | | 33 Hz | |
| Wavelength, nm | 1064/1053 nm | | | | 1053 nm | 1064/1053 nm | 1064 nm |
| Pulse energy ³⁾ | 8/10 mJ | 16/20 mJ | 32/40 mJ | 60/65 mJ | 70 mJ | 40 mJ | 60 mJ |
| Typical pulse duration ⁴⁾ | < 7 ns | | | < 5 ns | | | |
| Beam divergence ⁵⁾ | < 1.5 mrad | | | < 1 mrad | | | |
| Typical beam diameter ⁶⁾ | 1.5 mm | | 3.0 mm | 4.0 mm | | | |
| Pulse energy stability | < 0.5 % RMS | | | | | | |
| Power drift | ± 3.0 % | | | | | | |
| Beam profile | bell-shaped, > 80 % fit to Gaussian | | | | | | |
| Jitter | < 0.5 ns RMS | | | | | | |
| Polarization | linear, horizontal | | | | | | |

Optional harmonics generator ⁷⁾

| Pulse energy ³⁾ | | | | | | | |
|----------------------------|------------|----------|----------|----------|--------|-------|--------|
| 532/526.5 nm | 4/5 mJ | 8/10 mJ | 16/20 mJ | 30/32 mJ | 35 mJ | 20 mJ | 30 mJ |
| 355/351 nm | 2.4/3 mJ | 4.5/6 mJ | 9/12 mJ | 18/20 mJ | 21 mJ | 12 mJ | 18 mJ |
| 266/263 nm | 1.2/1.5 mJ | 2.4/3 mJ | 5/6 mJ | 9/10 mJ | 11 mJ | 6 mJ | 9 mJ |
| 213/211 nm | 0.5 mJ | 1 mJ | 2 mJ | 3.5 mJ | 4.2 mJ | 2 mJ | 2.5 mJ |



Optional attenuator ⁸⁾

| | | | | | | | |
|--------------------|--------|--|--|--|--|--|--|
| Transmission range | 1–95 % | | | | | | |
|--------------------|--------|--|--|--|--|--|--|

Operating requirements

| | | | | | | | |
|---------------------------|--|------|------|------|--|-------|--|
| Cooling requirements | air-cooled (water-free) | | | | | | |
| Ambient temperature | 15–30 °C | | | | | | |
| Relative humidity | 10–80 % (non-condensing) | | | | | | |
| Mains voltage | 90–230 VAC, single phase, 47–63 Hz ⁹⁾ | | | | | | |
| Average power consumption | 30 W | 40 W | 70 W | 80 W | | 100 W | |

- Due to continuous improvements all specifications are subject to change. Unless stated otherwise all specifications are measured at fundamental wavelength and maximum pulse repetition rate. The parameters marked typical are not specifications. They are indications of typical performance and will vary with each unit we manufacture.
- Factory-set pulse repetition rate is fixed at max repetition rate shown in the table. In internal triggering mode repetition rate can be divided by integer number down to f/2, f/3, f/4,... 1 Hz.
- When pulse energy is presented in xx/yy format, first number is for 1064 nm wavelength version, second – for 1053 nm version.
- At FWHM level at fundamental wavelength, measured with 350 ps rise time photodiode. Short pulse duration version is available, with pulse duration shorter by approx 50 %. Inquire for detailed specifications.
- Full angle measured at the 4 σ level.
- Beam diameter is measured 20 cm from laser output at the 4 σ level.
- Q2 is compatible with our attachable second harmonic generator (model SHG) and all models of stand-alone H-SMART harmonics generator. Pulse energies presented here are maximum values. Please refer to harmonic generator datasheets for detailed specifications.
- Motorized attenuator intended to be attached to the laser housing. Transmission can be changed remotely through laser web-server control interface.
- Laser can be powered from an appropriate 12 or 28 VDC power source. Please inquire for details.



Laser head Q2 and H-SMART harmonics generator
with FC fiber coupler on BSP2 base plate

SPECIFICATIONS at 50 – 200 Hz pulse repetition rate ¹⁾

| Model | Q2 | | | | | | | | |
|--------------------------------------|-------------------------------------|------------|----------|---------|----------|--------|------------|--------|----------|
| | -A50 | -B50 | -C50 | -D50 | -200 | -100 | -A100 | -B100 | -C100 |
| Pulse repetition rate ²⁾ | 50 Hz | | | | 200 Hz | 100 Hz | | | |
| Wavelength | 1064/1053 nm | | | 1064 nm | | | | | |
| Pulse energy ³⁾ | 5 mJ | 8/10 mJ | 16/20 mJ | 40 mJ | 1 mJ | 2.5 mJ | 5 mJ | 10 mJ | 20 mJ |
| Typical pulse duration ⁴⁾ | < 7 ns | | < 6 ns | < 5 ns | < 10 ns | | < 8 ns | | < 7 ns |
| Beam divergence ⁵⁾ | < 1 mrad | < 1.5 mrad | < 1 mrad | | < 2 mrad | | < 1.5 mrad | | < 1 mrad |
| Typical beam diameter ⁶⁾ | 1.5 mm | | 2.5 mm | 3.5 mm | 1.5 mm | | 2.0 mm | 2.5 mm | 3.5 mm |
| Pulse energy stability | < 0.5 % RMS | | | | | | | | |
| Power drift | ± 3.0 % | | | | | | | | |
| Beam profile | bell-shaped, > 80 % fit to Gaussian | | | | | | | | |
| Jitter | < 0.5 ns RMS | | | | | | | | |
| Polarization | linear, horizontal | | | | | | | | |

Optional harmonics generator ⁷⁾

| Pulse energy ³⁾ | | | | | | | | | |
|----------------------------|--------|------------|----------|----------|---------|---------|---------|--------|-------|
| 532/526.5 nm | 2.5 mJ | 4/5 mJ | 8/10 mJ | 16/20 mJ | 0.5 mJ | 1.25 mJ | 2.5 mJ | 5 mJ | 10 mJ |
| 355/351 nm | 1.5 mJ | 2.4/3 mJ | 4.5/6 mJ | 9/12 mJ | 0.25 mJ | 0.7 mJ | 1.5 mJ | 3 mJ | 6 mJ |
| 266/263 nm | 0.7 mJ | 1.2/1.5 mJ | 2.4/3 mJ | 5/6 mJ | 0.1 mJ | 0.3 mJ | 0.7 mJ | 1.5 mJ | 3 mJ |
| 213/211 nm | 0.2 mJ | 0.5 mJ | 1 mJ | 2 mJ | 0.02 mJ | 0.1 mJ | 0.25 mJ | 0.5 mJ | 1 mJ |

Optional attenuator ⁸⁾

| | |
|--------------------|--------|
| Transmission range | 1–95 % |
|--------------------|--------|


Operating requirements

| | |
|---------------------------|--|
| Cooling requirements | air-cooled (water-free) |
| Ambient temperature | 15–30 °C |
| Relative humidity | 10–80 % (non-condensing) |
| Mains voltage | 90–230 VAC, single phase, 47–63 Hz ⁹⁾ |
| Average power consumption | 30 W 50 W 80 W 100 W 40 W 50 W 70 W 80 W 100 W |

DIMENSIONS


| Unit | Dimensions (W×L×H) | Weight |
|-----------------------------|---------------------------------|--------|
| Laser head Q2 | 160 × 230 × 141 mm ³ | 4.8 kg |
| Controller unit | 108 × 172 × 59 mm ³ | 0.8 kg |
| Power adapter ¹⁾ | 50 × 125 × 31 mm ³ | 1.3 kg |
| Harmonic generator | 160 × 242 × 141 mm ³ | 4–6 kg |
| SHG module | 153 × 78 × 65 mm ³ | 0.7 kg |
| BSP2 plate | 219 × 491 × 32 mm ³ | 3 kg |

1. Power adapter dimensions might differ from indicated here, depending on model.




- Due to continuous improvements all specifications are subject to change. Unless stated otherwise all specifications are measured at fundamental wavelength and maximum pulse repetition rate. The parameters marked typical are not specifications. They are indications of typical performance and will vary with each unit we manufacture.
- Factory-set pulse repetition rate is fixed at max repetition rate shown in the table. In internal triggering mode repetition rate can be divided by integer number down to $f/2$, $f/3$, $f/4$,... 1 Hz.
- When pulse energy is presented in xx/yy format, first number is for 1064 nm wavelength version, second – for 1053 nm version.
- At FWHM level at fundamental wavelength, measured with 350 ps rise time photodiode. Short pulse duration version is available, with pulse duration shorter by approx 50 %. Inquire for detailed specifications.
- Full angle measured at the 4σ level.
- Beam diameter is measured 20 cm from laser output at the 4σ level.
- Q2 is compatible with our attachable second harmonic generator (model SHG) and all models of stand-alone H-SMART harmonics generator. Pulse energies presented here are maximum values. Please refer to harmonic generator datasheets for detailed specifications.
- Motorized attenuator intended to be attached to the laser housing. Transmission can be changed remotely through laser web-server control interface.
- Laser can be powered from an appropriate 12 or 28 VDC power source. Please inquire for details.


PERFORMANCE




Pulse duration of Q2-D10-1064



Beam profile of Q2-D10-1064
Axis length (X, Y) – (3.16, 3.18) mm,
effective diameter – 3.17 mm,
ellipticity – 99.2 %



Long term stability of Q2-D10-1064
Power drift: +0.1%; -0.1 %



Q2-F20-1053 1st harmonic beam profile
Axis length (X, Y) – (3.58, 3.47) mm,
effective diameter – 3.53 mm,
ellipticity – 96.9 %

PART NUMBERS

Q2-E20-1064-AT1

Model

Pulse energy level

- no letter → 0.1–4 mJ
- A → 5–7 mJ
- B → 8–15 mJ
- C → 16–30 mJ
- D → 31–50 mJ
- E → 51–69 mJ
- F → 70–90 mJ

Default pulse repetition rate in Hz

Optional items

- SHG0, SHG1, SHG2 → attachable second harmonic generator
- AT1 → motorized attenuator
- WH → water-cooled heatsink
- GL → low power CW guiding laser
- FC/NIR → fiber coupled output

Laser wavelength

- 1064 → Nd:YAG; 1064.2 nm
- 1053 → Nd:YLF; 1053.0 nm

OPTIONAL ITEMS

| | |
|--------------|---|
| WF | Stand-alone wireless router for wireless laser control |
| RS | Stand-alone adapter for laser control via RS-232 port |
| PC | Laptop computer for laser control |
| EXP | Stand-alone pulse generator for variable repetition rate |
| CT19 | 19" mounted controller with integrated AC/DC power supply |
| CTA19 | 19" mounted controller with air-purging unit |
| CTBR | Front/rear panel with brackets for standard controller |
| PS19 | 19" form factor AC/DC power supply |
| APU2 | Stand-alone air-purging unit with integrated AC/DC power supply |
| CST | Custom model |

IMAGES



Laser head Q2.
Front view



Laser head Q2.
Rear view



Fiber coupler FC



Motorized attenuator AT



Pulse energy monitor EM/BB



Laser controller CT



Laser controller CT.
Front view



Laser controller CT.
Rear view



Laser controller with air-purging unit CTA19




Front view of CTA19




Rear view of CTA19

DRAWINGS



Outline drawings of Q2 laser
(dimension in mm)



Outline drawings of Q2 laser
with AT motorized attenuator
(dimension in mm)